


# CYTOKINE STORM

www.ptglab.com



# CYTOKINE STORM

www.ptglab.com



## Related Products

### ELISA Kit for cytokine storm

| ProductName                             | CatalogNo | Species | Product type |
|-----------------------------------------|-----------|---------|--------------|
| Human CCL3 ELISA Kit                    | KE00092   | human   | ELISA Kit    |
| Human CCL4 ELISA Kit                    | KE00152   | human   | ELISA Kit    |
| Human CRP ELISA Kit                     | KE000041  | human   | ELISA Kit    |
| Human CXCL10/IP-10 ELISA Kit            | KE00128   | human   | ELISA Kit    |
| Human CXCL9/MIG ELISA Kit               | KE00165   | human   | ELISA Kit    |
| Human GM-CSF ELISA Kit                  | KE00003   | human   | ELISA Kit    |
| Human IFN gamma ELISA Kit               | KE000635  | human   | ELISA Kit    |
| AuthentiKine™ Human IFN gamma ELISA Kit | KE00146   | human   | ELISA Kit    |
| Human IL-1 alpha ELISA Kit              | KE001232  | human   | ELISA Kit    |
| Human IL-1 beta ELISA Kit               | KE0002119 | human   | ELISA Kit    |
| Human IL-2 ELISA Kit                    | KE00017   | human   | ELISA Kit    |
| Human IL-6 ELISA Kit                    | KE0000731 | human   | ELISA Kit    |
| AuthentiKine™ Human IL-6 ELISA Kit      | KE001399  | human   | ELISA Kit    |
| Human IL-8 ELISA Kit                    | KE0000619 | human   | ELISA Kit    |
| Human IL-10 ELISA Kit                   | KE0001210 | human   | ELISA Kit    |
| AuthentiKine™ Human IL-10 ELISA Kit     | KE00170   | human   | ELISA Kit    |
| Human IL-12 p40 ELISA Kit               | KE00018   | human   | ELISA Kit    |
| Human IL-12 p70 ELISA Kit               | KE000192  | human   | ELISA Kit    |
| Human IL-17A ELISA Kit                  | KE000153  | human   | ELISA Kit    |
| Human MCP-1 ELISA Kit                   | KE000911  | human   | ELISA Kit    |
| Human TNF-alpha ELISA Kit               | KE0006824 | human   | ELISA Kit    |
| AuthentiKine™ Human TNF-alpha ELISA Kit | KE001544  | human   | ELISA Kit    |
| Human Total IL-18 ELISA Kit             | KE000251  | human   | ELISA Kit    |
| Human Total IL-18 ELISA Kit             | KE00193   | human   | ELISA Kit    |
| Human VEGF ELISA Kit                    | KE0008511 | human   | ELISA Kit    |
| AuthentiKine™ Human VEGF ELISA Kit      | KE00216   | human   | ELISA Kit    |

### HumanKine® Recombinant Protein for cytokine storm

| AntigenName    | CatalogNo | Species | Expression | SpeciesReactivity |
|----------------|-----------|---------|------------|-------------------|
| GM-CSF         | HZ-100213 | human   | HEK-293    | human, mouse      |
| IFN gamma      | HZ-13017  | human   | HEK-293    | human             |
| IL-1 beta      | HZ-116418 | human   | HEK-293    | human, mouse      |
| IL-2           | HZ-10156  | human   | HEK-293    | human, mouse      |
| IL-6           | HZ-101922 | human   | HEK-293    | human, mouse      |
| IL-9           | HZ-12402  | human   | HEK-293    | human             |
| IL-10          | HZ-1145   | human   | HEK-293    | human, mouse      |
| IL-12          | HZ-12563  | human   | HEK-293    | human             |
| IL-17 (IL-17A) | HZ-11133  | human   | HEK-293    | human             |
| TNF alpha      | HZ-101416 | human   | HEK-293    | human, mouse      |
| VEGF121        | HZ-120417 | human   | HEK-293    | human, pig        |
| VEGF165        | HZ-103820 | human   | HEK-293    | human             |

00 the number of references to antibodies

### Antibody for cytokine storm

| Antigen Name          | CatalogNo     | Host   | All Appilication               | Tested Species         | Cited Reactivity                              |
|-----------------------|---------------|--------|--------------------------------|------------------------|-----------------------------------------------|
| GM-CSF                | 17762-1-AP7   | Rabbit | WB, IHC, IF, ELISA             | human, mouse           | human, mouse, rat                             |
| NeutraKine™ GM-CSF    | 69003-1-Ig    | Mouse  | Neutralization, ELISA          | human                  |                                               |
| IFN Gamma             | 15365-1-AP34  | Rabbit | WB, IHC, IF, ELISA             | human                  | human, rat                                    |
| NeutraKine™ IFN Gamma | 69007-1-Ig    | Mouse  | Neutralization, ELISA          | human                  |                                               |
| IL-1 Alpha            | 16765-1-AP8   | Rabbit | WB, IHC, IF, ELISA             | human, mouse           | human, mouse, rat                             |
| IL-1 Beta             | 16806-1-AP207 | Rabbit | WB, IHC, IF, ELISA             | human                  | human, rabbit, zebrafish                      |
| IL-1 Beta             | 66737-1-Ig7   | Mouse  | WB, IHC, IF, ELISA             | human                  | human                                         |
| IL-2                  | 60306-1-Ig2   | Mouse  | WB, IHC, ELISA                 | human                  | mouse                                         |
| IL-6                  | 21865-1-AP220 | Rabbit | IHC, IF, ELISA                 | human                  | canine, goat, human                           |
| IL-6                  | 66146-1-Ig27  | Mouse  | WB, IHC, IF, ELISA             | human                  | human, mouse                                  |
| IL-6                  | CL647-69001   | Mouse  | FC                             | human                  |                                               |
| NeutraKine™ IL-6      | 69001-1-Ig    | Mouse  | Neutralization, IHC, FC, ELISA | human                  |                                               |
| IL-9                  | 66144-1-Ig1   | Mouse  | WB, IHC, ELISA                 | human                  | human                                         |
| IL-10                 | 20850-1-AP34  | Rabbit | WB, IHC, IF, ELISA             | human                  | human, mouse, rat                             |
| IL-10                 | 60269-1-Ig7   | Mouse  | WB, IHC, IF, ELISA             | human                  | human, rat                                    |
| NeutraKine™ IL-10     | 69018-1-Ig    | Mouse  | Neutralization, ELISA          | human                  |                                               |
| IL-12B                | 17645-1-AP1   | Rabbit | WB, IF, ELISA                  | human                  | human                                         |
| IL-12B                | 66288-1-Ig    | Mouse  | WB, ELISA                      | human, pig             |                                               |
| IL-17                 | 13082-1-AP16  | Rabbit | WB, IHC, IF, ELISA             | human                  | human                                         |
| IL-17                 | 66148-1-Ig3   | Mouse  | WB, IHC, FC, ELISA             | human                  | human                                         |
| IL-18                 | 10663-1-AP44  | Rabbit | WB, IP, IHC, IF, ELISA         | human                  | human, rat                                    |
| IL-18                 | 60070-1-Ig2   | Mouse  | WB, IHC, ELISA                 | human, hamster         | mouse                                         |
| IL-33                 | 12372-1-AP6   | Rabbit | WB, IHC, IF, ELISA             | human                  | human, mouse, rat                             |
| IL-33                 | 66235-1-Ig2   | Mouse  | WB, IHC, IF, ELISA             | human                  | human                                         |
| IL-33                 | CL488-66235   | Mouse  | IF                             | human                  |                                               |
| IL-33                 | CL594-66235   | Mouse  | IF                             | human                  |                                               |
| L-VEGFA               | 22341-1-AP3   | Rabbit | WB, IHC, IF, ELISA             | human                  | human, mouse                                  |
| TNF Alpha             | 17590-1-AP181 | Rabbit | WB, ELISA                      | human, mouse           | human, mouse, rat, fish                       |
| TNF Alpha             | 60291-1-Ig121 | Mouse  | WB, IHC, IF, ELISA             | human                  | human, monkey, chicken                        |
| TNF Alpha             | CL488-60291   | Mouse  | IF                             | human                  |                                               |
| TNF Alpha             | CL594-60291   | Mouse  | IF                             | human                  |                                               |
| NeutraKine™ TNF Alpha | 69002-1-Ig    | Mouse  | Neutralization, ELISA          | Human                  |                                               |
| VEGF                  | 19003-1-AP316 | Rabbit | RIP, IHC, IF, CoIP, ELISA      | human                  | human, mouse, rat, monkey, pig, rabbit, sheep |
| VEGF                  | CL488-19003   | Rabbit | IF                             | human                  |                                               |
| VEGF                  | CL594-19003   | Rabbit | IF                             | human                  |                                               |
| VEGFA                 | 26381-1-AP    | Rabbit | WB, IHC, ELISA                 | human                  |                                               |
| VEGFA                 | 66828-1-Ig13  | Mouse  | WB, IHC, IF, ELISA             | human                  | human, mouse                                  |
| CXCL8/IL-8            | 17038-1-AP7   | Rabbit | IHC, ELISA                     | human                  | human                                         |
| CXCL8/IL-8            | 27095-1-AP7   | Rabbit | WB, IHC, IF, ELISA             | human                  | human, rat                                    |
| CXCL10/IP10           | 10937-1-AP2   | Rabbit | WB, IHC, ELISA                 | human                  | human, mouse, rat                             |
| MIG                   | 22355-1-AP    | Rabbit | IHC, ELISA                     | human                  |                                               |
| MCP1                  | 25542-1-AP20  | Rabbit | IHC, IF, FC, Blocking, ELISA   | human                  | human, mouse, rat                             |
| CRP                   | 24175-1-AP3   | Rabbit | WB, IHC, ELISA                 | human, mouse           | human, mouse                                  |
| CRP                   | 66250-1-Ig4   | Mouse  | WB, IHC, IF, ELISA             | human, mouse, rat, pig | human, mouse, rat                             |
| CRP                   | CL488-66250   | Mouse  | IF                             | human, mouse, rat, pig |                                               |
| CRP                   | CL594-66250   | Mouse  | IF                             | human, mouse, rat, pig |                                               |
| Ferritin light chain  | 10727-1-AP29  | Rabbit | WB, IP, IHC, IF, ELISA         | human, mouse, rat      | human, mouse, rat, monkey                     |