

For Research Use Only

SMAD4 Polyclonal antibody

Catalog Number: 10231-1-AP

Featured Product

102 Publications



Basic Information

Catalog Number:

10231-1-AP

Concentration:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0299

GenBank Accession Number:

BC002379

GeneID (NCBI):

4089

UNIPROT ID:

Q13485

Full Name:

SMAD family member 4

Calculated MW:

60 kDa

Observed MW:

63-70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

IP: 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC: 1:500-1:2000

Applications

Tested Applications:

WB, IHC, IP, ELISA

Cited Applications:

WB, IHC, IF, IP, CoIP, ChIP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat, pig, bovine, hamster, goat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: COLO 320 cells, HEK-293 cells, NIH/3T3 cells, mouse liver tissue, HCT 116 cells, HEK-293T cells, Neuro-2a cells

IP: HeLa cells,

IHC: human skin cancer tissue, human tonsillitis tissue, human cervical cancer tissue, mouse pancreas tissue, rat small intestine tissue

Background Information

Mammalian homologs of the Drosophila Mad gene include Smad1, Smad2, Smad3, Smad4 (DPC4), Smad5, Smad6, Smad7 and Smad8. Smad1 and Smad5 are effectors of BMP2 and BMP4 function while Smad2 and Smad3 are involved in TGF and activin-mediated growth modulation. Smad4 has been shown to mediate all of the above activities through interaction with various Smad family members. Smad6 and Smad7 regulate the response to activin/ TGF signaling by interfering with TGF-mediated phosphorylation of other Smad family members. This antibody is a rabbit polyclonal antibody raised against an internal region of human SMAD4.

Notable Publications

| Author | Pubmed ID | Journal | Application |
|----------|-----------|----------------|-------------|
| Lu Liu | 34576190 | Int J Mol Sci | WB, IP |
| Bing Li | 34528447 | Epigenomics | WB |
| Yanli Xu | 30206204 | Cell Death Dis | CoIP |

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 10231-1-AP (SMAD4 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.




IP result of anti-SMAD4 (IP:10231-1-AP, 4ug; Detection:10231-1-AP 1:1500) with HeLa cells lysate 1320 ug.




Immunohistochemical analysis of paraffin-embedded human skin cancer tissue slide using 10231-1-AP (SMAD4 antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded mouse pancreas tissue slide using 10231-1-AP (SMAD4 antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using 10231-1-AP (SMAD4 antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded rat small intestine tissue slide using 10231-1-AP (SMAD4 antibody) at dilution of 1:1000 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).