

For Research Use Only

MOG Polyclonal antibody

Catalog Number: 12690-1-AP

26 Publications



Basic Information

Catalog Number: 12690-1-AP	GenBank Accession Number: BC035938	Purification Method: Antigen affinity purification
Concentration: 500 ug/ml	GeneID (NCBI): 4340	Recommended Dilutions: WB: 1:500-1:3000
Source: Rabbit	UNIPROT ID: Q16653	IHC: 1:50-1:500
Isotype: IgG	Full Name: myelin oligodendrocyte glycoprotein	IF-P: 1:50-1:500
Immunogen Catalog Number: AG3273	Calculated MW: 295 aa, 34 kDa	
	Observed MW: 25-28 kDa	

Applications

Tested Applications: WB, IHC, IF-P, ELISA	Positive Controls:
Cited Applications: WB, IHC, IF	WB : mouse brain tissue, rat brain tissue
Species Specificity: human, mouse, rat	IHC : mouse brain tissue,
Cited Species: human, mouse, rat	IF-P : mouse brain tissue,
Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0	

Background Information

Myelin/oligodendrocyte glycoprotein (MOG), a 23-28 kDa glycoprotein, a myelin antigen at the outer surface of the central nervous system (CNS) myelin sheath, which may trigger T-cell as well as B-cell responses. It therefore constitutes a pivotal target for autoimmune responses, which result in inflammation and also demyelination in the CNS. Its presence on the outer-most lamellae of mature CNS myelin and its late appearance during myelinogenesis suggest that it contributes to myelin maturation or maintenance. 10 isoforms of MOG produced by alternative splicing have been described, and heterodimers may be formed between the different isoforms. Defects in MOG are the cause of narcolepsy type 7 (NRCLP7), a neurological disabling sleep disorder characterized by excessive daytime sleepiness, sleep fragmentation, symptoms of abnormal rapid-eye-movement (REM) sleep, cataplexy, hypnagogic hallucinations, and sleep paralysis. Role of MOG in the pathogenesis of multiple sclerosis (MS) has been reported but remains to be clarified.

Notable Publications

Author	Pubmed ID	Journal	Application
Isabella Farhy-Tselnick	34494546	Elife	IF
Alessandro Dinoto	36257153	Mult Scler Relat Disord	WB
Simona Perga	33051914	Brain Pathol	IHC

Storage

Storage: Store at -20°C. Stable for one year after shipment.
Storage Buffer: PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 12690-1-AP (MOG antibody) at dilution of 1:1500 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 12690-1-AP (MOG antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunofluorescent analysis of (4% PFA) fixed paraffin-embedded mouse brain tissue using MOG antibody (12690-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).