

For Research Use Only

CNOT7 Polyclonal antibody

Catalog Number: 14102-1-AP

Featured Product

9 Publications



Basic Information

Catalog Number:

14102-1-AP

Concentration:

253 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5235

GenBank Accession Number:

BC060852

GeneID (NCBI):

29883

UNIPROT ID:

Q9UIV1

Full Name:

CCR4-NOT transcription complex,
subunit 7

Calculated MW:

285 aa, 33 kDa

Observed MW:

33-35 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:1000

IP: 0.5-4.0 µg for 1.0-3.0 mg of total
protein lysate

Applications

Tested Applications:

WB, IP, ELISA

Cited Applications:

WB, IP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB: mouse testis tissue,

IP: mouse testis tissue,

Background Information

Cnot7, Ccr4-Not complex 7, is a component of CCR4-NOT complex and known to be a transcription factor and a modulator of mRNA degradation in yeast. Mammalian Cnot7 interacts with Tob and its family members, which is known as an inhibitor of BMP signaling through interfering with the Smad system and is involved in bone metabolism [PMID:17451368]. It is one of at least 8 subunits that make up the CCR4-Not transcription complex, which is involved in negative regulation of the pol II directed transcription in yeast²⁴ and may be involved in global gene regulation by modulation of TATAA-binding protein (TBP) activity [PMID:10880468].

Notable Publications

Author	Pubmed ID	Journal	Application
Masashi Yamaji	28297718	Nature	WB
Hsin-Hui Hsieh	34090459	J Inflamm (Lond)	WB
Roy Matkovic	35013187	Nat Commun	IP

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 14102-1-AP (CNOT7 antibody) at dilution of 1:500 incubated at room temperature for 1.5 hours.



IP result of anti-CNOT7 (IP:14102-1-AP, 4ug; Detection:14102-1-AP 1:300) with mouse testis tissue lysate 4800ug.