

For Research Use Only

ASNS Polyclonal antibody

Catalog Number: 14681-1-AP

Featured Product

61 Publications



Basic Information

Catalog Number:

14681-1-AP

Concentration:

600 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6393

GenBank Accession Number:

BC014621

GeneID (NCBI):

440

UNIPROT ID:

P08243

Full Name:

asparagine synthetase

Calculated MW:

64 kDa

Observed MW:

55-64 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

IP: 0.5-4.0 µg for 1.0-3.0 mg of total protein lysate

IHC: 1:50-1:500

IF/ICC: 1:50-1:500

Applications

Tested Applications:

IHC, IF/ICC, IP, ELISA

Cited Applications:

WB, IHC, IF

Species Specificity:

human

Cited Species:

human, mouse, rat

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

IP: HeLa cells,

IHC: human retinoblastoma tissue, human gliomas tissue

IF/ICC: K-562 cells,

Background Information

Asparagine synthetase [glutamine-hydrolyzing](ASNS) is also named as TS11. The asparagine synthetase gene, encoding the enzyme that catalyzes the synthesis of asparagine and glutamate using glutamine and aspartate, is a gene for which transcription is highly regulated by the nutritional status of the cell(PMID:15385533). ASNS has some isoforms with the MW of 55-64 kDa.

Notable Publications

Author	Pubmed ID	Journal	Application
César Vásquez-Trincado	34550363	Hum Mol Genet	WB
Shiqi Wang	34515223	Pharmacology	WB
Jie Jiang	31659118	J Biol Chem	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




WB result of ASNS antibody (14681-1-AP; 1:6000; incubated at room temperature for 1.5 hours) with sh-Control and sh-ASNS transfected Jurkat cells.




IP result of anti-ASNS (IP:14681-1-AP, 4ug; Detection:14681-1-AP 1:500) with HeLa cells lysate 1200ug.




Immunohistochemical analysis of paraffin-embedded human retinoblastoma tissue slide using 14681-1-AP (ASNS Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human retinoblastoma tissue slide using 14681-1-AP (ASNS Antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Various lysates were subjected to SDS PAGE followed by western blot with 14681-1-AP (ASNS antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (4% PFA) fixed K-562 cells using ASNS antibody (14681-1-AP) at dilution of 1:200 and Multi-rAb CoraLite® Plus 488-Goat Anti-Rabbit Recombinant Secondary Antibody (H+L) (RGAR002).