

For Research Use Only

HNRNPA2B1 Polyclonal antibody

Catalog Number: 14813-1-AP

Featured Product

71 Publications



Basic Information

Catalog Number:

14813-1-AP

Concentration:

400 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6563

GenBank Accession Number:

BC000506

GeneID (NCBI):

3181

UNIPROT ID:

P22626

Full Name:

heterogeneous nuclear
ribonucleoprotein A2/B1

Calculated MW:

37 kDa

Observed MW:

36 kDa and 38 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:5000-1:50000

IHC: 1:1000-1:4000

Applications

Tested Applications:

WB, IHC, ELISA

Cited Applications:

WB, IHC, IF, IP, CoIP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat, pig, chicken

**Note-IHC: suggested antigen retrieval with
TE buffer pH 9.0; (*) Alternatively, antigen
retrieval may be performed with citrate
buffer pH 6.0**

Positive Controls:

WB : A549 cells, A431 cells, HL-60 cells, HSC-T6 cells,
HeLa cells, K-562 cells

IHC : human colon tissue, human lymphoma tissue,
human colon cancer tissue

Background Information

HnRNP A2B1 is one of the heterogeneous nuclear ribonucleoproteins (HnRNPs), which regulate the splicing and transportation of mRNA. It has a role in packaging nascent mRNA, alternative splicing, cytoplasmic RNA trafficking, translation, and stabilization. Besides it may involve in telomere maintenance, cell proliferation and differentiation, and glucose transport. Gene HURNPA2B1 encodes two proteins, HNRNPA2 (36kd) and HNRNPB1 (38kd), by alternative splicing.

Notable Publications

Author	Pubmed ID	Journal	Application
Huimin Guo	33134163	Front Oncol	WB
Feiya Zhu	34631545	Front Oncol	IHC
Xin Liu	34496888	Mol Cancer	WB,IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



WB result of HNRNPA2B1 antibody (14813-1-AP; 1:3000; incubated at room temperature for 1.5 hours) with sh-Control and sh-HNRNPA2B1 transfected A549 cells.



Various lysates were subjected to SDS PAGE followed by western blot with 14813-1-AP (HNRNPA2B1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human normal colon slide using 14813-1-AP (HNRNPA2B1 antibody) at dilution of 1:2000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).