

For Research Use Only

# DNMT3L Polyclonal antibody

Catalog Number: 14939-1-AP

2 Publications



## Basic Information

Catalog Number:

14939-1-AP

Size:

300 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6822

GenBank Accession Number:

BC002560

GeneID (NCBI):

29947

UNIPROT ID:

Q9UJW3

Full Name:

DNA (cytosine-5-)-methyltransferase

3-like

Calculated MW:

44 kDa

Observed MW:

44 kDa and 70 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:1000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse

Cited Species:

human

Positive Controls:

WB : mouse testis tissue, PC-3 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
Lijie Zhou	34185414	Clin Transl Med	WB
Yang Fu	36508844	J Reprod Immunol	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse testis tissue were subjected to SDS PAGE followed by western blot with 14939-1-AP (DNMT3L Antibody) at dilution of 1:600 incubated at room temperature for 1.5 hours.