

For Research Use Only

OAS1/3 Polyclonal antibody

Catalog Number:14955-1-AP

Featured Product

33 Publications



Basic Information

Catalog Number:

14955-1-AP

Concentration:

500 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6793

GenBank Accession Number:

BC000562

GeneID (NCBI):

4938

UNIPROT ID:

P00973

Full Name:

2',5'-oligoadenylate synthetase 1,
40/46kDa

Calculated MW:

46 kDa

Observed MW:

~40 kDa, 100-120 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

IHC: 1:100-1:400

IF/ICC: 1:200-1:800

Applications

Tested Applications:

WB, IHC, IF/ICC, ELISA

Cited Applications:

WB, IHC, IF, CoIP

Species Specificity:

human, mouse

Cited Species:

human, mouse, rat, pig

**Note-IHC: suggested antigen retrieval with
TE buffer pH 9.0; (*) Alternatively, antigen
retrieval may be performed with citrate
buffer pH 6.0**

Positive Controls:

WB: IFN alpha treated HeLa cells, mouse brain tissue

IHC: human liver tissue,

IF/ICC: A431 cells,

Background Information

The 2-prime,5-prime oligoadenylate synthetases (OASs), such as OAS1, are interferon-induced proteins characterized by their capacity to catalyze the synthesis of 2-prime,5-prime oligomers of adenosine (2-5As). OAS1 (type I enzymes) has some isoforms with the MW of 40, 42, 44, 46, and 48 kDa (PMID:12590567, 19383565). OAS1 is a strong candidate for determining susceptibility or resistance to viral infections(PMID:15732009). This antibody may have cross-reaction with OAS3 protein due to the high homology.

Notable Publications

Author	Pubmed ID	Journal	Application
Huanru Wang	31575039	Int J Mol Sci	WB
Arthur Wickenhagen	34581622	Science	WB
Li Che	34731560	Acta Biochim Pol	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Untreated and IFN alpha treated HeLa cells were subjected to SDS PAGE followed by western blot with 14955-1-AP (OAS1/3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 14955-1-AP (OAS1/3 Antibody) at dilution of 1:200 (under 10x lens).



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using 14955-1-AP (OAS1/3 Antibody) at dilution of 1:200 (under 40x lens).



Immunofluorescent analysis of (-20°C Methanol) fixed A431 cells using OAS1/3 antibody (14955-1-AP) at dilution of 1:400 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).