

For Research Use Only

# GABARAP Polyclonal antibody

Catalog Number: 18723-1-AP

Featured Product

28 Publications



## Basic Information

### Catalog Number:

18723-1-AP

### Concentration:

900 ug/ml

### Source:

Rabbit

### Isotype:

IgG

### GenBank Accession Number:

BC106748

### GeneID (NCBI):

11337

### UNIPROT ID:

O95166

### Full Name:

GABA(A) receptor-associated protein

### Calculated MW:

14 kDa

### Observed MW:

14-18 kDa

### Purification Method:

Antigen affinity purification

### Recommended Dilutions:

WB: 1:1000-1:6000

IP: 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IF/ICC: 1:50-1:500

FC (Intra): 0.40 ug per 10<sup>6</sup> cells in a 100 µl suspension

## Applications

### Tested Applications:

WB, IF/ICC, FC (Intra), IP, ELISA

### Cited Applications:

WB, IF, CoIP

### Species Specificity:

human, mouse

### Cited Species:

human, mouse, monkey

### Positive Controls:

**WB** : HeLa cells, A549 cells, mouse brain tissue, human placenta tissue, mouse pancreas tissue, human heart tissue, mouse skeletal muscle tissue, HuH-7 cells, SH-SY5Y cells, RAW 264.7 cells, mouse liver tissue, rat brain tissue

**IP** : HeLa cells,

**IF/ICC** : HepG2 cells, Starvation treated HepG2 cells

**FC (Intra)** : HepG2 cells,

## Background Information

GABARAP (Gamma-aminobutyric acid type A receptor-associated protein), was initially identified as GABAA receptor binding protein and has shown to be implicated in intracellular trafficking of GABAA receptors. Through interacting with a broad range of proteins, GABARAP has been demonstrated as a multifunctional protein implicated in intracellular vesicle trafficking, apoptosis, and autophagy. The protein level of lipidated form of GABARAP (GABARAP-II) is upregulated during the autophagic process, which makes it as an autophagic marker. This antibody detected endogenous GABARAP protein in lysates from various tissues. It may cross-react with GABARAPL1 due to the high identity of amino acid sequence.

## Notable Publications

Author	Pubmed ID	Journal	Application
Linyu Sun	3455354	Mol Cell	WB
Svenja Zielke	33111629	Autophagy	WB
Zhen-Hua Chen	27740626	Cell Death Differ	WB

## Storage

### Storage:

Store at -20°C. Stable for one year after shipment.

### Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 18723-1-AP (GABARAP antibody) at dilution of 1:3000 incubated at room temperature for 1.5 hours.



1X10<sup>6</sup> HepG2 cells were intracellularly stained with 0.4 ug Anti-Human GABARAP (18723-1-AP) and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) at dilution 1:1000 (red), or 0.4 ug Isotype Control. Cells were fixed with 4% PFA and permeabilized with Flow Cytometry Perm Buffer (PF00011-C).



Immunofluorescent analysis of (-20°C Ethanol) fixed Chloroquine treated HepG2 cells and untreated HepG2 cells using GABARAP antibody (18723-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).



IP result of anti-GABARAP (IP:18723-1-AP, 4ug; Detection:18723-1-AP 1:2000) with HeLa cells lysate 1200 ug.