

For Research Use Only

# TMEM41A Polyclonal antibody

Catalog Number: 20768-1-AP

Featured Product

3 Publications



## Basic Information

Catalog Number:

20768-1-AP

Size:

400 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG14753

GenBank Accession Number:

BC008043

GeneID (NCBI):

90407

UNIPROT ID:

Q96HV5

Full Name:

transmembrane protein 41A

Calculated MW:

264 aa, 30 kDa

Observed MW:

26-28 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

## Applications

Tested Applications:

WB, ELISA

Cited Applications:

WB

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse

Positive Controls:

WB : mouse liver tissue, mouse small intestine tissue, rat liver tissue, HepG2 cells

## Background Information

### Notable Publications

Author	Pubmed ID	Journal	Application
H-Heinrich Hoffmann	33052348	bioRxiv	WB
H-Heinrich Hoffmann	33338421	Cell	WB
Keigo Morita	30093494	J Cell Biol	WB

## Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 20768-1-AP (TMEM41A antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.