

For Research Use Only

FAM134B Polyclonal antibody

Catalog Number: 21537-1-AP

Featured Product

51 Publications



Basic Information

Catalog Number:

21537-1-AP

Concentration:

800 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG16176

GenBank Accession Number:

BC030517

GeneID (NCBI):

54463

UNIPROT ID:

Q9H6L5

Full Name:

family with sequence similarity 134, member B

Calculated MW:

497 aa, 55 kDa

Observed MW:

65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:2000-1:16000

IP: 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC: 1:50-1:500

Applications

Tested Applications:

WB, IHC, IP, ELISA

Cited Applications:

WB, IHC, IF, IP

Species Specificity:

human, mouse, rat

Cited Species:

human, mouse, rat, pig, chicken

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB : C2C12 cells, HEK-293 cells, rat brain tissue, mouse brain, mouse cerebellum

IP : Jurkat cells,

IHC : mouse brain tissue, mouse heart tissue

Background Information

FAM134B, also known as the reticulophagy regulator 1 (RETREG1) or JK-1, is a member of the family with sequence similarity 134. FAM134B is an important ER-phagy receptor. FAM134B regulates the size and shape of the ER and participates in some ER-phagy-related processes (PMID: 33199694). FAM134B inhibition contributes to impair proteostasis in the ER due to the accumulation of misfolded or aggregated proteins, which in turn leads to compromised neuronal survival and progressive neuronal degenerative diseases. In addition, FAM134B mutations are common in patients with colorectal adenocarcinoma, and oesophageal squamous cell carcinoma (PMID: 29226326).

Notable Publications

Author	Pubmed ID	Journal	Application
Mark L Schultz	30202070	Nat Commun	WB
Smriti Parashar	34467852	Elife	WB, IF
Svenja Zielke	33111629	Autophagy	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




WB result of FAM134B antibody (21537-1-AP; 1:2000; incubated at room temperature for 1.5 hours) with sh-Control and sh-FAM134B transfected HEK-293 cells.



Various lysates were subjected to SDS PAGE followed by western blot with 21537-1-AP (FAM134B antibody) at dilution of 1:8000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded mouse brain tissue slide using 21537-1-AP (FAM134B antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



IP result of anti-FAM134B (IP:21537-1-AP, 4ug; Detection:21537-1-AP 1:5000) with Jurkat cells lysate 1600 ug.