

For Research Use Only

LY6E Polyclonal antibody, PBS Only

Catalog Number: 22144-1-PBS

Featured Product



Basic Information

Catalog Number:

22144-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG17573

GenBank Accession Number:

BC119708

GeneID (NCBI):

4061

UNIPROT ID:

Q16553

Full Name:

lymphocyte antigen 6 complex, locus E

Calculated MW:

131 aa, 14 kDa

Observed MW:

12-16 kDa

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Background Information

The predicted MW of LY6E is 13 kDa. Catalog#22144-1-AP recognises two band of 13 kDa and 18 kDa, and the additional 18 kDa band may be due to glycosylation.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 22144-1-AP (LY6E antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 22144-1-PBS in a different storage buffer formulation.