

For Research Use Only

FOXA2 Polyclonal antibody

Catalog Number: 22474-1-AP

Featured Product

25 Publications



Basic Information

Catalog Number:

22474-1-AP

Concentration:

550 ug/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18221

GenBank Accession Number:

BC011780

GeneID (NCBI):

3170

UNIPROT ID:

Q9Y261

Full Name:

forkhead box A2

Calculated MW:

457 aa, 48 kDa

Observed MW:

55-65 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB: 1:500-1:2000

IP: 0.5-4.0 ug for 1.0-3.0 mg of total protein lysate

IHC: 1:2000-1:8000

IF-P: 1:50-1:500

Applications

Tested Applications:

WB, IHC, IF-P, IP, ELISA

Cited Applications:

WB, IHC, IF, CoIP, ChIP

Species Specificity:

human

Cited Species:

human, mouse

Positive Controls:

WB: HT-29 cells, HepG2 cells, A549 cells, MDA-MB-231 cells

IP: HepG2 cells,

IHC: human colon cancer tissue,

IF-P: rat stomach tissue,

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Background Information

FOXA2 also named as HNF 3 beta or TCF3B is a 457 amino acid protein, which contains one fork-head DNA-binding domains. FOXA2 can Shuttle between the nucleus and cytoplasm in a CRM1-dependent manner. Phosphorylation on FOXA2 Thr-156 abolishes binding to target promoters and subsequent transcription activation upon INS stimulation. FOXA2 as a transcription factor involves in embryonic development, establishment of tissue-specific gene expression and regulation of gene expression in differentiated tissues. The molecular weight of FOXA2 is 48kDa, but the phosphorylated FOXA2 may be a 55-65kDa protein.

Notable Publications

Author	Pubmed ID	Journal	Application
Lin Yang	31754329	Int J Biol Sci	WB
Yan Huang	36334352	Biomaterials	WB,IF
Sang Luo	34926672	Ann Transl Med	IF

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data




Various lysates were subjected to SDS PAGE followed by western blot with 22474-1-AP (FOXA2 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.




Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 22474-1-AP (FOXA2 antibody) at dilution of 1:4000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using 22474-1-AP (FOXA2 antibody) at dilution of 1:4000 (under 20x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).



IP result of anti-FOXA2 (IP:22474-1-AP, 4ug; Detection:22474-1-AP 1:500) with HepG2 cells lysate 1320 ug.



Immunofluorescent analysis of (4% PFA) fixed rat stomach tissue using FOXA2 antibody (22474-1-AP) at dilution of 1:200 and CoraLite®488-Conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).