

For Research Use Only

SARS-COV-2 NSP1 Polyclonal antibody



Catalog Number: 29287-1-AP

Basic Information

Catalog Number:

29287-1-AP

Size:

1000 μ g/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG30935

GenBank Accession Number:

YP_009742608

GeneID (NCBI):

43740578

Full Name:

ORF1ab

Purification Method:

Antigen affinity purification

Applications

Tested Applications:

ELISA

Species Specificity:

virus

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SARS-CoV-2 Envelope protein Antibody (29287-1-AP) tested by ELISA. SARS-COV-2 NSP1 protein was coated onto microtiter plates at 1 µg/well and then incubated with a dilution series of SARS-COV-2 NSP1 protein Antibody (29287-1-AP). Bound antibodies were detected with HRP conjugated anti-Rabbit IgG followed by incubation with HRP Substrate and then measuring the resulting absorbance at 450 nm.