

For Research Use Only

# COXIV Monoclonal antibody

Catalog Number: 60251-1-Ig **10 Publications**



## Basic Information

Catalog Number: 60251-1-Ig	GenBank Accession Number: BC021236	Purification Method: Protein G purification
Concentration: 1000 µg/ml	GeneID (NCBI): 1327	CloneNo.: 3C7D2
Source: Mouse	UNIPROT ID: P13073	
Isotype: IgG1	Full Name: cytochrome c oxidase subunit IV isoform 1	
Immunogen Catalog Number: AG20551	Calculated MW: 19.6 kDa	
	Observed MW: 17-18 kDa	

## Applications

Tested Applications:  
ELISA

Cited Applications:  
WB, IF

Species Specificity:  
human

Cited Species:  
human, mouse, rat

## Background Information

COX4l1, also named as COX4 and COXIV-1, belongs to the cytochrome c oxidase IV family. It is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport. COX4l1 is a marker for mitochondria. It has two isoforms (isoform 1 and 2). Isoform 1(COX4l1) is ubiquitously expressed and isoform 2 is highly expressed in lung tissues. COX4l1 is commonly used as a loading control. This antibody is specific to COX4l1 and do not cross reacts with COX4l2.

## Notable Publications

Author	Pubmed ID	Journal	Application
Zhijiang Chen	35726572	Eur J Histochem	WB,IF
Bo Zhang	36588386	Autophagy	WB
Yufei Song	29459731	Sci Rep	IF

## Storage

Storage:  
Store at -20°C. Stable for one year after shipment.

Storage Buffer:  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:  
T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

