

For Research Use Only

APOA5 Monoclonal antibody

Catalog Number: 65015-1-Ig



Basic Information

Catalog Number:

65015-1-Ig

Size:

500 µg/ml

Source:

Mouse

Isotype:

IgG1

GenBank Accession Number:

BC101787

GeneID (NCBI):

116519

UNIPROT ID:

Q6Q788

Full Name:

apolipoprotein A-V

Calculated MW:

366 aa, 41 kDa

Observed MW:

39 kDa-41 kDa

Purification Method:

Antigen affinity purification

Recommended Dilutions:

WB 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : Human serum,

Background Information

Storage

Storage:

Store at 2-8°C.

Storage Buffer:

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

