

For Research Use Only

Aequorea victoria GFP tag Monoclonal antibody, PBS Only

Catalog Number: 66002-1-PBS



Basic Information

Catalog Number:

66002-1-PBS

Concentration:

1 mg/ml

Source:

Mouse

Isotype:

IgG2a

Immunogen Catalog Number:

AG2128

GenBank Accession Number:

M62653

GeneID (NCBI):

Full Name:

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

1E10H7

Applications

Tested Applications:

WB, IF/ICC, IP, ELISA

Species Specificity:

recombinant protein

Background Information

Green Fluorescent Proteins (GFPs) encompass a diverse range of proteins carrying a green chromophore, originating from various species and forming different protein lineages. Wildtype GFP consists of 238 amino acid residues (26.9 kDa). GFP was first identified in the jellyfish *Aequorea victoria*. It emits green light with a peak wavelength of 509 nm upon excitation by blue light at 395 nm. When fused with other proteins, GFP serves as a versatile reporter protein e.g. for quantifying expression levels or facilitates visualization of subcellular localization through fluorescence microscopy. This antibody is a mouse (IgG2a) monoclonal antibody against GFP reactive to GFP, eGFP, eYFP, YFP and CFP.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



HEK293 cells, CHO cells, normal mouse brain and GFP transgenic mouse brain were subjected to SDS PAGE followed by western blot with 66002-1-Ig (GFP tag antibody) at dilution of 1:200000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 66002-1-PBS in a different storage buffer formulation.



IP result of anti-GFP tag (IP:66002-1-Ig, 4ug; Detection:66002-1-Ig 1:10000) with Transfected HEK-293T cells lysate 400 ug. This data was developed using the same antibody clone with 66002-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (-20°C Ethanol) fixed Transfected HEK-293 cells using GFP tag antibody (66002-1-Ig, Clone: 1E10H7) at dilution of 1:800 and CoraLite® 594-Conjugated AffiniPure Goat Anti-Mouse IgG(H+L) (SA00013-3). This data was developed using the same antibody clone with 66002-1-PBS in a different storage buffer formulation.