

AIM2 Monoclonal antibody

Catalog Number: 66902-1-Ig

Featured Product

18 Publications

Basic Information

Catalog Number:

66902-1-Ig

Concentration:

2000 ug/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG14607

GenBank Accession Number:

BC010940

GeneID (NCBI):

9447

UNIPROT ID:

O14862

Full Name:

absent in melanoma 2

Calculated MW:

343 aa, 39 kDa

Observed MW:

37-39 and 53 kDa

Purification Method:

Protein G purification

CloneNo.:

3E12A8

Recommended Dilutions:

WB: 1:1000-1:4000

IHC: 1:500-1:2000

Applications

Tested Applications:

WB, IHC, ELISA

Cited Applications:

WB, IHC, IF, IP, CoIP

Species Specificity:

human, mouse, rat, rabbit

Cited Species:

human, mouse

Note-IHC: suggested antigen retrieval with TE buffer pH 9.0; (*) Alternatively, antigen retrieval may be performed with citrate buffer pH 6.0

Positive Controls:

WB: RAW 264.7 cells, rabbit spleen tissue, human heart tissue, Sp2/O cells

IHC: human urothelial carcinoma tissue, rat small intestine tissue

Background Information

Absent in melanoma 2 (AIM2) is a sensor of cytosolic DNA that is responsible for activation of the inflammasome and host immune responses to DNA viruses and intracellular bacteria [PMID:22695634]. Upon binding to DNA, it is thought to undergo oligomerization and to associate with PYCARD initiating the recruitment of caspase-1 precursor and processing of interleukin-1 beta and interleukin-18. Otherwise, it functions as a tumor suppressor which may act by repressing NF-kappa-B transcriptional activity [PMID:16432157,20566831]. 66902-1-Ig antibody detects the 37-39 and 53 kDa protein in SDS-PAGE similar to the paper published. [PMID: 28809949]

Notable Publications

Author	Pubmed ID	Journal	Application
Ruby E Dawson	34489308	Gut	CoIP
Hao Liu	36245295	EMBO J	WB,IF
Qingyu Liang	34815793	Theranostics	WB

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage


For technical support and original validation data for this product please contact:

T: 4006900926


E: Proteintech-CN@ptglab.comW: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 66902-1-Ig (AIM2 antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



Immunohistochemical analysis of paraffin-embedded human urothelial carcinoma tissue slide using 66902-1-Ig (AIM2 antibody) at dilution of 1:1000 (under 10x lens). Heat mediated antigen retrieval with Tris-EDTA buffer (pH 9.0).