

FOR IN VITRO RESEARCH USE ONLY.
NOT FOR USE IN HUMANS OR ANIMALS.

Recombinant human HSP90AB1 protein



Basic Information

Catalog Number:

Ag1947

Size:

50 µg

Form:

Available lyophilized

Species:

human

Expression Source:

e coli-derived, PGEX-4T, with N-terminal GST.

Biological Activity:

Not tested

Endotoxin Level:

Please contact the lab for more information

Peptide Sequence:

RALLFIPRRAPFDLFENKKKKNNIKLYVRRVFIMDSCDE
LIPEYLNFI RGVVDS E D L P L N I S R E M L Q Q S K I L K V I R K N I
V K K C L E L F S E L A E D K E N Y K K F Y E A F S K N L K L G I H E D S T N
R R R L S E L L R Y H T S Q S G D E M T S L S E Y V S R M K E T Q K S I Y Y I
T G E S K E Q V A N S A F V E R V R K R G F E V V Y M T E P I D E Y C V
Q Q L K E F D G K S L V S V T K E G L E P E D E E E K K K M E E S K A K F
E N L C K L M K E I L D K K V E K V T I S N R L V S S P C C I V T S T Y G W
T A N M E R I M K A Q A L R D N S T M G Y M M A K K H L E I N P D H P I
V E T L R Q K A E A D K N D K A V K D L V L L F E T A L L S S G F S L E D
P Q T H S N R I Y R M I K L G L G I D E D E V A A E E P N A A V P D E I P P
(330-710 aa encoded by BC014485)

Selected Validation Data

For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.