

For Research Use Only

CoraLite®594-conjugated TLR8 Monoclonal antibody

Catalog Number:CL594-67317

1 Publications



Basic Information

Catalog Number:

CL594-67317

Concentration:

1000 µg/ml

Source:

Mouse

Isotype:

IgG1

Immunogen Catalog Number:

AG16566

GenBank Accession Number:

BC101075

GeneID (NCBI):

51311

UNIPROT ID:

Q9NR97

Full Name:

toll-like receptor 8

Calculated MW:

1041 aa, 120 kDa

Observed MW:

120 kDa

Purification Method:

Protein G purification

CloneNo.:

3B6A4

Recommended Dilutions:

IF-P: 1:50-1:500

Excitation/Emission maxima wavelengths:

588 nm / 604 nm

Applications

Tested Applications:

IF-P

Cited Applications:

IF

Species Specificity:

Human

Cited Species:

mouse

Positive Controls:

IF-P: human tonsillitis tissue,

Background Information

Toll-like receptors (TLRs) are a group of pathogen recognition receptors (PRRs) that can be divided into two subfamilies depending on their subcellular localization(PMID: 32268138). Toll-like receptor 8 (TLR8), a sensor for pathogen-derived single-stranded RNA (ssRNA), binds to uridine (Uri) and ssRNA to induce defense responses (PMID:30535046). TLR8 proteins were expressed in spleen, epididymal sperm, testicular sperm cells, and Leydig cells but were weak in Sertoli cells. It can identify two isoforms,120 and 122 kDa.

Notable Publications

Author	Pubmed ID	Journal	Application
Ga-Yeon Son	39479800	Function (Oxf)	IF

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite®594 TLR8 antibody (CL594-67317, Clone: 3B6A4) at dilution of 1:200.



Immunofluorescent analysis of (4% PFA) fixed human tonsillitis tissue using CoraLite®594 TLR8 antibody (CL594-67317, Clone: 3B6A4) at dilution of 1:200.