

For Research Use Only.
Not For Use In Diagnostics.

Streptavidin-CoraLite® Plus 700 Conjugate

Catalog Number: CL700-PF00030



Information

Catalog Number:
CL700-PF00030
Host:
Applications:
FC, IF

Reactivity:
Physical State:
Liquid
Conjugation:
CoraLite®700

Recommended Dilutions

1:200-1:1000 for flow cytometry and immunofluorescence

Storage

Storage:
2-8°C for one year.
Storage Buffer:
PBS with 0.1% sodium azide, pH 7.4.
Aliquoting is unnecessary for -20°C storage


Purity

For technical support and original validation data for this product please contact:

T: 4006900926 E: Proteintech-CN@ptglab.com W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Flow cytometric analysis of human PBMCs using Biotin anti-human CD8 and CoraLite®700-conjugated Streptavidin. Human PBMCs were stained with Biotin CD8 Recombinant Antibody (Cat.No. Biotin-65559, Clone: UCHT4), washed with PBS, and stained with CoraLite®700-conjugated Streptavidin (Cat.No.CL700-PF00030) and analyzed by flow cytometry.