

For Research Use Only

FITC Plus Anti-Mouse CD44 (IM7)

Catalog Number: FITC-65117

4 Publications



Basic Information

Catalog Number:

FITC-65117

Concentration:

100ug, 0.5 mg/mL

Source:

Rat

Isotype:

IgG2b, kappa

GenBank Accession Number:

BC051388

GeneID (NCBI):

12505

UNIPROT ID:

P15379

Full Name:

CD44 antigen

Purification Method:

Affinity purification

CloneNo.:

IM7

Recommended Dilutions:

FC: 0.5 ug per 10⁶ cells in 100 μ l suspension

Excitation/Emission maxima wavelengths:

495 nm / 524 nm

Applications

Tested Applications:

FC

Cited Applications:

FC

Species Specificity:

Mouse

Cited Species:

human, mouse

Positive Controls:

FC : mouse splenocytes,

Background Information

CD44 is a type I transmembrane glycoprotein expressed on embryonic stem cells and in various levels on other cell types including connective tissues and bone marrow. CD44 expression is also upregulated in subpopulations of cancer cells and is recognized as a molecular marker for cancer stem cells (PMID: 29747682). It is a cell-surface receptor that mediates cell-cell and cell-matrix interactions through its affinity for hyaluronic acid (HA) and possibly also through its affinity for other ligands (PMID: 10694938). Adhesion with HA plays an important role in cell migration, tumor growth and progression. CD44 is also involved in lymphocyte activation, recirculation and homing, and in hematopoiesis.

Notable Publications

Author	Pubmed ID	Journal	Application
Kunxiang Gong	39680126	Cell Mol Life Sci	FC
Yao Qi	39603539	J Control Release	FC
Kunxiang Gong	38733440	Discov Oncol	FC

Storage

Storage:

Store at 2-8°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 0.09% sodium azide, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



1X10⁶ mouse splenocytes were surface co-stained with Anti-Mouse CD3 and 0.5 ug FITC Anti-Mouse CD44 (FITC-65117, Clone:IM7) or 0.5 ug Control Antibody. Cells were not fixed.