

For Research Use Only

HRP-conjugated GFP tag Monoclonal antibody

Catalog Number:HRP-66002

19 Publications



Basic Information

Catalog Number:

HRP-66002

Concentration:

1000 ug/ml

Source:

Mouse

Isotype:

IgG2a

GenBank Accession Number:

M62653

GeneID (NCBI):

UNIPROT ID:

P42212

Full Name:

Calculated MW:

26 kDa

Observed MW:

-

Purification Method:

Antigen affinity purification

CloneNo.:

1E10H7

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB

Cited Applications:

WB, IF, CoIP

Species Specificity:

aequorea victoria, recombinant protein

Positive Controls:

WB : Recombinant protein,

Background Information

Green Fluorescent Proteins (GFPs) encompass a diverse range of proteins carrying a green chromophore, originating from various species and forming different protein lineages. Wildtype GFP consists of 238 amino acid residues (26.9 kDa). GFP was first identified in the jellyfish *Aequorea victoria*. It emits green light with a peak wavelength of 509 nm upon excitation by blue light at 395 nm. When fused with other proteins, GFP serves as a versatile reporter protein e.g. for quantifying expression levels or facilitates visualization of subcellular localization through fluorescence microscopy. This antibody is a HRP conjugated mouse (IgG2a) monoclonal antibody against GFP reactive to GFP, eGFP, eYFP, YFP and CFP. Note: Do not add Azium (Sodium Azide or Smite) into the dilution buffer. Azium is the HRP inhibitor which decreases the enzymatic activity of HRP.

Notable Publications

Author	Pubmed ID	Journal	Application
Shu-Yi Yang	31659125	Plant Physiol	
Zhihui Ruan	36274707	Front Microbiol	WB, CoIP
Xingyu Zhang	34018922	Elife	WB

Storage

Storage:

Store at -20°C. Avoid exposure to light. Stable for one year after shipment.

Storage Buffer:

PBS with 50% glycerol, 0.05% Proclin300, 0.5% BSA, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Recombinant protein were subjected to SDS PAGE followed by western blot with HRP-66002 (GFP tag Antibody) at dilution of 1:10000.