

IHCeasy[®] CLTA Ready-To-Use IHC Kit

Catalog Number: **KHC0592**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

Clathrin is the major protein of the polyhedral coat of coated pits and vesicles which entrap specific macromolecules during receptor-mediated endocytosis. The clathrin molecule has a triskelion shape. Each clathrin triskelion is composed of three identical heavy chains and three light chains of two types, LCA (CLTA) and LCB (CLTB). The light chain subunits are thought to regulate the formation or disassembly of clathrin coats.


Synonyms

Clathrin light chain A, clathrin, light chain (Lca), CLTA, LCA


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using KHC0592 (CLTA IHC Kit).




Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using KHC0592 (CLTA IHC Kit).



Immunohistochemical analysis of paraffin-embedded human kidney tissue slide using KHC0592 (CLTA IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using KHC0592 (CLTA IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat liver tissue slide using KHC0592 (CLTA IHC Kit).