

IHCeasy[®] RHOA Ready-To-Use IHC Kit

Catalog Number: **KHC0638**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

RhoA is a member of the Rho family of small GTPases, which cycle between inactive GDP-bound and active GTP-bound states and function as molecular switches in signal transduction cascades. Rho proteins promote reorganization of the actin cytoskeleton and regulate cell shape, attachment, and motility. Overexpression of RhoA is associated with tumor cell proliferation and metastasis. RhoA signalling is critical to many cellular processes including migration, mechanotransduction, and is often disrupted in carcinogenesis.


Synonyms

ARH12, ARHA, Rho cDNA clone 12, RHO12, RHOA, RHOH12, Transforming protein RhoA


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC0638 (RHOA IHC Kit).




Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using KHC0638 (RHOA IHC Kit).



Immunohistochemical analysis of paraffin-embedded human liver tissue slide using KHC0638 (RHOA IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse stomach tissue slide using KHC0638 (RHOA IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat stomach tissue slide using KHC0638 (RHOA IHC Kit).