

IHCeasy[®] G6PD Ready-To-Use IHC Kit

Catalog Number: **KHC0876**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

G6PD(Glucose-6-phosphate 1-dehydrogenase) belongs to the glucose-6-phosphate dehydrogenase family. It is a housekeeping enzyme encoded in mammals by an X-linked gene. G6PD has important functions in intermediary metabolism because it catalyzes the first step in the pentose phosphate pathway and provides reductive potential in the form of NADPH.


Synonyms

G6P, G6PD, G6PD1


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human lung cancer tissue slide using KHC0876 (G6PD IHC Kit).




Immunohistochemical analysis of paraffin-embedded human stomach cancer tissue slide using KHC0876 (G6PD IHC Kit).




Immunohistochemical analysis of paraffin-embedded human ovary tumor tissue slide using KHC0876 (G6PD IHC Kit).




Immunohistochemical analysis of paraffin-embedded human liver cancer tissue slide using KHC0876 (G6PD IHC Kit).



Immunohistochemical analysis of paraffin-embedded human colon cancer tissue slide using KHC0876 (G6PD IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using KHC0876 (G6PD IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using KHC0876 (G6PD IHC Kit).