

IHCeasy[®] CIRBP Ready-To-Use IHC Kit

Catalog Number: **KHC1218**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

CIRBP, also named as A18HNRNP and CIRP, is a cold-inducible mRNA binding protein that plays a protective role in the genotoxic stress response by stabilizing transcripts of genes involved in cell survival. CIRBP is also involved in cap-independent translation upon moderate cold-shock. It acts as a translational activator. CIRBP is a suppressive rather stimulatory effect on proliferation. CIRBP can translocate from nucleus to cytoplasm under stresses such as UV irradiation or heat shock. Suppression of CIRBP increased the apoptotic cell population of neural stem cells at moderate low temperature.


Synonyms

A18 hnRNP, A18HNRNP, CIRBP, CIRP


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human thyroid cancer tissue slide using KHC1218 (CIRBP IHC Kit).



Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using KHC1218 (CIRBP IHC Kit).




Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using KHC1218 (CIRBP IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse pancreas tissue slide using KHC1218 (CIRBP IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using KHC1218 (CIRBP IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat pancreas tissue slide using KHC1218 (CIRBP IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat colon tissue slide using KHC1218 (CIRBP IHC Kit).