

IHCeasy[®] CDC25A Ready-To-Use IHC Kit

Catalog Number: **KHC2228**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	


Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

CDC25A belongs to the MPI phosphatase family. It is a tyrosine protein phosphatase which functions as a dosage-dependent inducer of mitotic progression. CDC25A directly dephosphorylates CDC2 and stimulates its kinase activity. It also dephosphorylates CDK2 in complex with cyclin E, in vitro. CDC25A overexpression in human cancers might contribute to tumorigenesis.


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human cervical cancer tissue slide using KHC2228 (CDC25A IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse testis tissue slide using KHC2228 (CDC25A IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse colon tissue slide using KHC2228 (CDC25A IHC Kit).



Immunohistochemical analysis of paraffin-embedded mouse thymus tissue slide using KHC2228 (CDC25A IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat colon tissue slide using KHC2228 (CDC25A IHC Kit).