

IHCeasy[®] HDAC4 Ready-To-Use IHC Kit

Catalog Number: **KHC2330**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Rabbit Polyclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Rabbit

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

HDAC4, also named as HDACA, belongs to the histone deacetylase family. HDAC4 is responsible for the deacetylation of lysine residues on the N-terminal part of the core histones (H2A, H2B, H3 and H4). Histone deacetylation gives a tag for epigenetic repression and plays an important role in transcriptional regulation, cell cycle progression and developmental events. Histone deacetylases act via the formation of large multiprotein complexes. HDAC4 is involved in muscle maturation via its interaction with the myocyte enhancer factors such as MEF2A, MEF2C and MEF2D.


Synonyms

KIAA0288, HD4, EC:3.5.1.98


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2330 (HDAC4 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse small intestine tissue slide using KHC2330 (HDAC4 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse lung tissue slide using KHC2330 (HDAC4 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse salivary gland tissue slide using KHC2330 (HDAC4 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat lung tissue slide using KHC2330 (HDAC4 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat small intestine tissue slide using KHC2330 (HDAC4 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat salivary gland tissue slide using KHC2330 (HDAC4 IHC Kit).