

IHCeasy[®] OGDH Ready-To-Use IHC Kit

Catalog Number: **KHC2386**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

OGDH(2-oxoglutarate dehydrogenase, mitochondrial) is also named as OGDC-E1 and belongs to the alpha-ketoglutarate dehydrogenase family. It is one of the three enzymatic components of 2-oxoglutarate dehydrogenase complex, which catalyzes the overall conversion of 2-oxoglutarate to succinyl-CoA and CO₂.


Synonyms

EC:1.2.4.2, E1 α , Alpha-KGDH-E1, Alpha-ketoglutarate dehydrogenase, 2-oxoglutarate dehydrogenase, mitochondrial


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human placenta tissue slide using KHC2386 (OGDH IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse heart tissue slide using KHC2386 (OGDH IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat adrenal gland tissue slide using KHC2386 (OGDH IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat cerebellum tissue slide using KHC2386 (OGDH IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat heart tissue slide using KHC2386 (OGDH IHC Kit).