

IHCeasy[®] NLGN1 Ready-To-Use IHC Kit

Catalog Number: **KHC2670**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

NLGN1, also named KIAA1070, belongs to the type-B carboxylesterase/lipase family. Neuroligin1 (NLGN1) is a main component of the excitatory glutamatergic synapses complex and plays a role in synapse assembly and function.


Synonyms

NLG1,NLGN1,KIAA1070,neuroligin 1,Neuroligin-1


Selected Validation Data



Immunohistochemical analysis of paraffin-embedded mouse adrenal gland tissue slide using KHC2670 (NLGN1 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat adrenal gland tissue slide using KHC2670 (NLGN1 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat brain tissue slide using KHC2670 (NLGN1 IHC Kit).