

IHCeasy[®] SMCR7L/MID51 Ready-To-Use IHC Kit

Catalog Number: **KHC2942**

General Information

Sample type:
FFPE tissue
Cited sample type:
Reactivity:
Human, Mouse, Rat
Cited Reactivity:

Assay type:
Immunohistochemistry
Primary antibody type:
Mouse Monoclonal
Secondary antibody type:
Polymer-HRP-Goat anti-Mouse

Kit Component

Component	Size	Concentration
Antigen Retrieval Buffer	100 mL	50×
Washing Buffer	100 mL ×2	20×
Blocking Buffer	5 mL	RTU
Primary Antibody	5 mL	RTU
Secondary Antibody	5 mL	RTU
Chromogen Component A	0.2 mL	RTU
Chromogen Component B	4 mL	RTU
Signal Enhancer	5 mL	RTU
Counter Staining Reagent	5 mL	RTU
Mounting Media	5 mL	RTU
Control Slide	1 slide (Optional)	FFPE
Datasheet	1 Copy	
Manual	1 Copy	

Storage Instructions

All the reagents are stored at 2-8°C. The kit is stable for 6 months from the date of receipt.


Background

Human SMCR7L gene encodes, MID51, the mitochondrial dynamic protein of 51 kDa (also called mitochondrial elongation factor 1, MIEF1). MID51 is a single-pass membrane protein anchored to the mitochondrial outer membrane and regulates mitochondrial morphology. Mitochondrial morphology is controlled by two opposing processes: fusion and fission. Elevated MID51 levels induce extensive mitochondrial fusion, whereas depletion of MID51 causes mitochondrial fragmentation. MID51 interacts with and recruits Drp1 to mitochondria, suggesting a critical role of MID51 in regulation of mitochondrial fusion-fission machinery in vertebrates.


Synonyms

SMCR7L, MID51, MIEF1, SMCR7-like protein, Smith-Magenis syndrome chromosomal region candidate gene 7 protein-like


Selected Validation Data




Immunohistochemical analysis of paraffin-embedded human pancreas cancer tissue slide using KHC2942 (SMCR7L/MID51 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse kidney tissue slide using KHC2942 (SMCR7L/MID51 IHC Kit).




Immunohistochemical analysis of paraffin-embedded mouse liver tissue slide using KHC2942 (SMCR7L/MID51 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat heart tissue slide using KHC2942 (SMCR7L/MID51 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat kidney tissue slide using KHC2942 (SMCR7L/MID51 IHC Kit).



Immunohistochemical analysis of paraffin-embedded rat liver tissue slide using KHC2942 (SMCR7L/MID51 IHC Kit).