

For Research Use Only.  
Not For Use In Diagnostics.

# Streptavidin-PE Conjugate

Catalog Number: PE-PF00030



[www.ptglab.com](http://www.ptglab.com)

---

## Information

Catalog Number:

PE-PF00030

Host:

Applications:

FC, IF

Reactivity:

Physical State:

Liquid

Conjugation:

PE

---

## Recommended Dilutions

1:200-1:1000 for flow cytometry and immunofluorescence

---

## Storage

Storage:

2-8°C for one year.

Storage Buffer:

PBS with 0.1% sodium azide and 0.5% BSA, pH 7.4.

Aliquoting is unnecessary for -20°C storage

---

## Purity

---

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com) W: [ptgcn.com](http://ptgcn.com)


---

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Flow cytometric analysis of human PBMCs using Biotin anti-human CD3 and PE Streptavidin. Human PBMCs were stained with Biotin CD3 Recombinant Antibody (Cat.No. Biotin-65569, Clone: OKT3), washed with PBS, and stained with PE Streptavidin (Cat.No. PE-PF00030) and analyzed by flow cytometry.



Immunofluorescent analysis of (4% PFA) fixed HepG2 cells using Biotin PEX14 antibody (Biotin-10594) at dilution of 1:200 and PE Streptavidin (Cat.No. PE-PF00030), CL488-phalloidin (green).