

For Research Use Only

CoraLite®Plus 405链霉亲和素偶联物

Catalog Number: PF00028



产品介绍

该CoraLite®Plus 405链霉亲和素偶联物（CoraLite®Plus 405-conjugated Streptavidin）为溶液状态，通过流式细胞分析进行了质量控制。对于流式细胞染色，建议每次测试使用量为0.025 ug/test，务必每次实验都进行浓度优化以确保得到最佳染色。

产品成分

名称	CoraLite®Plus 405链霉亲和素偶联物
浓度	0.5 mg/ml
存储溶液	PBS with 0.1% sodium azide, pH 7.4.
激发波长	399 nm
发射波长	422 nm

包装规格

100 ug/500 ug/1 mg


保存条件

2°C-8°C避光保存，一年有效。

注意

CoraLite®Plus 405链霉亲和素偶联物无需稀释保存，2°C-8°C避光储存，不要冷冻。

Validation Data



Flow cytometric analysis of human PBMCs using Biotin anti-human CD3 and CoraLite® Plus 405-conjugated Streptavidin: human PBMCs were stained with Biotin CD8 Recombinant Antibody (Cat.No. Biotin-65559, Clone: UCHT4), washed with PBS, and stained with CoraLite® Plus 405-conjugated Streptavidin (Cat.No. PF00028) and... analyzed by flow cytometry.

For technical support and original validation data for this product please contact

T: 027-87531629

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.