

For Research Use Only

GFP-Trap® Magnetic Agarose, Kit for Immunoprecipitation



www.ptgcn.com

Catalog Number: gtmak

50 Publications

Basic Information

Catalog Number:
gtmak

Applications:
IP, CoIP, ChIP, RIP

Conjugate:
Magnetic agarose beads; bead size: ~40 µm (cross-linked 6 % magnetic agarose beads)

Host:
Alpaca

Type:
Nanobody

Class:
Recombinant

Description

The ChromoTek GFP-Trap® Magnetic Agarose are affinity beads for IP of GFP-fusion proteins. It comprises a GFP Nanobody/VHH coupled to magnetic agarose beads.

Binding capacity

15-20 µg of recombinant GFP per 25 µL bead slurry

Specificity/Target

AcGFP, Clover, eGFP, Emerald, GFP, GFP5, GFP Envy, GFP S65T, mGFP, mPhluorin, PA-GFP, Superfolder GFP, TagGFP, TagGFP2, monomeric eGFP A206K, CFP, YFP, Citrine, eCitrine, eYFP, Venus, Ypet, BFP
For the complete list, please click here: [Fluorescent protein specificity table](#)

Elution buffer

SDS sample buffer
0.2 M glycine pH 2.5

Affinity (K_D)

Dissociation constant K_D of 1 pM

Storage

Storage:
Shipped at ambient temperature. Upon receipt store at 4°C. Stable for one year. Do not freeze!
Storage Buffer:
20% ethanol


For technical support and original validation data for this product please contact:

T: 4006900926

E: Proteintech-CN@ptglab.com
W: www.ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



GFP-Trap Magnetic Agarose for immunoprecipitation of GFP-fusion proteins. Note single-band in bound fraction, which indicates a very effective pull-down. I: Input, FT: Flow-Through, B: Bound.