

For Research Use Only

# MK2-Trap Agarose, Kit for Immunoprecipitation



[www.ptgcn.com](http://www.ptgcn.com)

Catalog Number: **mtak**

## Basic Information

**Catalog Number:**

mtak

**Applications:**

IP, CoIP

**Conjugate:**

Agarose beads; bead size: ~ 90 µm (cross-linked 4 % agarose beads)

**Host:**

Alpaca

**Type:**

Nanobody

**Class:**

Recombinant

## Description

ChromoTek's MK2-Trap Agarose kit are affinity beads for IP of MK2-fusion proteins. It comprises a MK2 Nanobody/ VHH conjugated to agarose beads.

## Specificity/Target

MK2 (MAP kinase-activated protein kinase 2, MAPKAP-K2)

## Elution buffer

SDS sample buffer  
0.2 M glycine pH 2.5

## Storage

**Storage:**

Upon receipt store at +4°C. Do not freeze!

**Storage Buffer:**

20% ethanol

For technical support and original validation data for this product please contact:

T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [www.ptgcn.com](http://www.ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data

